

Cat. No.	A1L	A2L	A3L	A4L
Code	801L-012	801L-014	801L-018	801L-022
ISO No.	001-012	001-014	001-018	001-022
Head Size	012	014	018	022
Total Length	24	24	24	24
Fine				
Medium	▽	▽	▽	▽
Coarse		▽	▽	

**Surgical Length Round**

Code	Size
A1LM	Medium 012
A2LM	Medium 014
A2LC	Coarse 014
A3LM	Medium 018
A3LC	Coarse 018
A4LM	Medium 022

6 Pieces per Pack




Cat. No.	A2XL	A3XL
Code	801XL-014	801XL-018
ISO No.	001-014	001-018
Head Size	014	018
Total Length	28	28
Fine		
Medium	▽	▽
Coarse		

**Extra Long Surgical**

Code	Size
A2XLM	Medium 014
A3XLM	Medium 018


2 Pieces per Pack



Cat. No.	A4PF	A4SF	A4XXF
Code	801-022	801-022	801-022
ISO No.	001-022	001-022	001-022
Head Size	022	022	022
Total Length	19	19	19
Extra Extra Fine			▽
Super Fine		▽	
Magic Touch	▽		

**Round**

Code	Size
A4PF	Magic Touch 022 2/Pk
A4SF	Super Fine 022 6/Pk
A4XXF	Extra Extra Fine 022 3/Pk




Cat. No.	E3R
Code	889-008
ISO No.	169-008
Head Size	008
Head Length	3
Total Length	24
S. Fine	
Fine	▽
Medium	▽
Coarse	
S. Coarse	

**Separation**

Code	Size
E3RF	Fine 008
E3RM	Medium 008

6 Pieces per Pack



Cat. No.	F1RXL
Code	850-016
ISO No.	199-016
Head Size	016
Head Length	10
Total Length	31
S. Fine	
Fine	▽
Medium	▽
Coarse	▽
S. Coarse	

**Round End Taper**


Code	Size
F1RXLF	Fine 016
F1RXLM	Medium 016

2 Pieces per Pack



SEE DEMO VIDEO


The F1RXL is an extra-long 31mm bur, which offers fantastic clear visibility, allowing quicker, far more precise work. Designed by Dr. Amir Weizmann.



Cat. No.	PR2SF	PR2XXF
Code	390-018	390-018
ISO No.	254-018	254-018
Head Size	018	018
Total Length	19	19
Extra Extra Fine		▽
Super Fine	▽	

**PEAR**

Code	Size
PR2SF	Super Fine 018 6/Pk
PR2XXF	Extra Extra Fine 018 3/Pk




Cat. No.	PR1210	PR1910
Code	802L-016	802L-019
ISO No.	464-016	464-019
Head Size	012/016	018/019
Head Length	10	10
Total Length	23	23
S. Fine		
Fine		
Medium		
Coarse	▽	▽
S. Coarse		

**Endo Access**

Code	Size
PR1210C	Coarse 012/016
PR1910C	Coarse 018/019

2 Pieces per Pack

A special endodontic bur, combining a round bur for initial penetration with a Diamond plated long taper to widen the pulp chamber and the canal walls



Cat. No.	L4
Code	851-014
ISO No.	220-014
Head Size	014
Head Length	9
Total Length	24
S. Fine	
Fine	
Medium	▽
Coarse	
S. Coarse	

**Safe End**

Code	Size
L4M	Medium 016

2 Pieces per Pack



Cat. No.	M3
Code	856L-018
ISO No.	198-018
Head Size	018
Head Length	9
Total Length	21
S. Fine	
Fine	
Medium	▽
Coarse	
S. Coarse	

**Round End Taper**

Code	Size
M3M	Medium 018
C1C	Coarse 014
C0C	Coarse 012

6 Pieces per Pack

## Minimal Invasive Burs

Minimal invasive Diamond burs for opening of fissures / endodontic use.



Cat. No.	MA0M	MA1M	MA2M	MA25M
Code	801-009	801-012	801-014	801-016
ISO No.	001-009	001-012	001-014	001-016
Head Size	9	12	14	16
Head Length				
Total Length	21	22	22	22
Grit Size	▽	▽	▽	▽

Code	Size
MA0M	Medium 009
MA1M	Medium 012
MA2M	MEDIUM 014
MA25M	Medium 016
3 Pieces per Pack	



Cat. No.	MAB0M	MAB3M
Code	830-008	830-010
ISO No.	237-008	237-010
Head Size	8	10
Head Length	2.5	2.5
Total Length	19	19
Grit Size	▽	▽

Code	Size
MAB0M	Medium 008
MAB3M	Medium 010
3 Pieces per Pack	



### ENDODONTIC KITS

### COMPLETE ENDODONTIC BUR KIT



ORDER No. ENDOSET1

Profile	Cat. No.	Application
	A4PF	Safe Ceramic penetration without microfractures
	FGB331	Easy & fast Metal penetration
	FG-ENDOZ152	Pulp chamber opening without perforation risk
	E3RM	Metal post removal / Metal screw weakening cementation
	E3RF	Safe search of canal entrance / span
	A1LM	Finding the canal entrance / span
	A4LM	Pulp chamber opening
	PR2SF	Middling Composite polishing
	PR2XXF	Gentle Composite polishing & finishing
	A4SF	Middling Composite polishing
	A4XXF	Gentle Composite polishing & finishing




# Endo Z

## Endo Access

Endo Safe End Bur

- Specially designed Carbide bur with non-cutting safe-end tip.
- Safely widens pulp chamber without perforation and ledging.
- Creates funnel shape to allow easier access to root canals.
- Safe ended tip prevents penetration of the pulp chamber floor.



Cat. No.	EndoZ
Head Size	016
Head Length	9
Total Length	23
FG	✳

Code	Size
ENDOZ	016 Non cutting Tip
3 Pieces per Pack	

## ||||| Barracuda Metal Cutting

## ||||| Carbide Bur



Cat. No.	B1558
Head Size	012
Head Length	4.1
FG	✳

Round End Cross Cut Fissure

Code	Size
B1558	012
3 Pieces per Pack	



Cat. No.	B330	B331
Head Size	008	010
Head Length	2	2
FG	✳	✳


Pear

Code	Size
FGB330	008
FGB331	010
3 Pieces per Pack	



Cat. No.	4
Head Size	014
Head Length	1.1
FG	
FGSS	
FGOS	✳

Code	Size
FG4OS	014
3 Pieces per Pack	



### Munce Discovery Burs® Troughers

**The Shallow Troughers (31mm)** were specifically designed for procedures on the chamber floor or just beyond such as troughing for the MB2 canal of maxillary molars, troughing for the middle mesial canal of mandibular molars or to “brush out” the isthmus between the MB & ML canals of mandibular molars. Shallow Troughers are also useful where limited inter-arch space precludes the use of the Deep Troughers.

**The Deep Troughers (34mm)** were designed for deeper operations within the root such as exposure of separated instruments, removal of pulpal tissue from deep within C-shaped canals and chasing severely-calcified canals into the mid-root and beyond.

**The Cariesectomy Bur (31mm)** has a #6 round carbide head and was designed to facilitate the type of deep caries removal often required in endodontic cases, and for which there was previously no comparable instrument. This bur—often in combination with the #4 Shallow Trougher—is ideal for caries removal from around crown margins from within the access cavity.

**Features**

- Slow-speed
- Stiff Shaft
- No Heat Generation
- Maintains a Functional View Corridor beyond the Handpiece Head
- Latch-type
- Durable Carbide Head

Code	Size
Munce Discovery Burs® 31mm Shallow Troughers	
S0VP	Variety pack 6 Pieces/Pk
S4025	# ¼ Pink
S4050	# ½ Grey
S4100	# 1 Purple
S4200	# 2 White
S4300	# 3 Yellow
S4400	# 4 Red
S4600	# 6 Green Cariesectomy Bur
4 Pieces per Pack	

Code	Size
Munce Discovery Burs® 34mm Deep Troughers	
D0VP	Variety pack 6 Pieces/Pk
D4025	# ¼ Pink
D4050	# ½ Grey
D4100	# 1 Purple
D4200	# 2 White
D4300	# 3 Yellow
D4400	# 4 Red
4 Pieces per Pack	