


# TOMIDENT PRODUCT CATALOGUE 2024



#ThinkEndoThinkTomident

# 2024 PRODUCT CATALOGUE





**| 6**    **PATENTED  
ROTARY FILES**

**| 8**    **ROTARY FILES**

**| 12**   **RECIPROCATION  
FILES**

**| 14**   **GLIDE PATH  
FILES**

**| 14**   **SPECIAL  
ROTARY FILES**

**| 16**   **HAND FILES**

**| 18**   **SEALERS,  
LINERS & TEMP  
MATERIALS**

**| 20**   **GUTTA PERCHA  
& PAPER POINTS**

**| 21**   **LENTULOS,  
GATES GLIDDEN  
& PEESO REAMER**

**| 22**   **MTA**

**| 24**   **CORE BUILDUP  
& POSTS**

**| 25**   **IRRIGATION  
SYRINGES , NEEDLES  
& RUBBER DAM**

**| 26**   **EQUIPMENT**

**| 38**   **BURS**

**| 41**   **INSTRUMENTS**

**| 42**   **SCALER TIPS**

# Welcome to Tomident.

## Our Story

An Australian family owned company focused on providing new, innovative and cost effective products and solutions.

We are committed to be the industry leader through continued support and education.



Tomident was started with the aim of bringing quality service and customer experience values back to our industry.

We are dedicated to focus on you the customer and listening to your feedback about the products you demand from around the world.

With over 100 years combined industry experience we believe our team is able to effectively support you and your team.

A handwritten signature in black ink, appearing to read 'Nicholas Varkas', with a long horizontal line extending to the right.

Nicholas Varkas

Managing Director

All Trademarks, Registered Names, or Brands belong to their respective owners. Tomident does not intend to mislead or sell replica, copy or remanufacture goods. Comparable product references are intended to guide you to a better, similar or complementary item.

All items supplied are TGA registered unless otherwise stated with the phrase "Sim use".

The information appearing in our website, promotional material, email or catalogues are of general nature and is NOT a substitute for, or to be used in conjunction with professional medical advice. The products we supply are for Dental Professionals with AHPRA registration.

Tomident reserved the right to change, update, cancel or substitute items, pricing or promotions at any time.

# TIME FOR A BETTER FILE ?

**BUY 10  
GET 10%  
OFF**  
MIX n MATCH  
ALL ROTARY

**Freedom to  
Choose Size  
& Tapers**


**What do I want  
from my new File  
System ?**

**Easy to  
follow  
Sequence**

**Add a  
Glider  
to your  
system**

**Medium  
Flexibility**

**AF BLUE ROTARY**



ORDER NO. A04... \$55  
PACK OF 6

Medium Flexibility

Triangular Cross Section

#17/12 Open File  
#18/05 Glide Path

#20/04 to #40/04  
Plus  
#20/06 to #40/06

**Highest  
Flexibility**

**F ONE**



ORDER NO. A20... \$55  
PACK OF 3+1

Two active cutting points


A Flexible file for all scenarios.

Can be used in a simple single file or multiple file sequence for more complexed Endo cases.

Pack of 3 Shaping Files + 1 Open file

**Reciprocating**

**AF BLUE R3**



ORDER NO. A05... \$55  
PACK OF 3


Popular Reciprocating 30° 150° movement

R3 Suited to a clinic that completes fewer or less complicated cases by choice.

Simplicity and predictability in One file

**Reciprocating**

**W+**



ORDER NO. W... \$66  
PACK OF 6


Popular Reciprocating 30° 150° movement

W+ Suited to a clinic that wants to continue using this popular file system but wants improved flexibility with a cost effective price.

STERILE BLISTER

**Glide Path**

**C-PATH FILE**



ORDER NO. F21... \$55  
PACK OF 6

Glide Path Sizes


#13/02  
#16/02  
#19/02

-Inactive Tip

Shape Memory Wire Technology

**Medium  
Flexibility**

**AF BLUE S4**



ORDER NO. A04... \$49.5  
PACK OF 4

Medium Flexibility Triangular Cross Section


Easy to follow 4 file sequence kit.

#17/12 Open File  
#18/05 Glide Path  
#25/06 Shaping File  
#35/04 Shaping File

Ability to expand system further with AF Blue Rotary

**High  
Flexibility**

**VTAPER GOLD**



ORDER NO. V01... V02... \$55  
PACK OF 6

Popular SX S1 S2 F1 F2 F3 sequence kit.

High Flexibility

Convex Triangular cross section

No need to change your existing workflow

**Highest  
Flexibility**

**TX+**



ORDER NO. TX... \$66  
PACK OF 6

Popular TXA PRG TX1 TX2 TX3 TX4 sequence kit.

Three Wire Technology Archimedean spiral cross section

No need to change your existing workflow

STERILE BLISTER

**Highest  
Flexibility**

**F ONE SIE KIT**




ORDER NO. SIE... \$55  
PACK OF 4


A Highly Flexible file in a Sequence kit for many scenarios.

#13/03 Glide Path  
#20/04 Shaping File  
#25/06 Shaping File  
#35/04 Shaping File

Ability to expand further with AF F One for complex cases.



Instructions for use  
AF FONE<sup>®</sup> FILE



Optional Files:  
Glide Path: 16/02, 19/02  
Working Length: 20/06, 25/06, 35/06



# PATENTED CROSS SECTION

## F ONE™

FLEXIBLE - EFFICIENT - FREEDOM OF CHOICE  
SIMPLE EASY TO FOLLOW SEQUENCES

- Patented Cross section
- Unique Silver Heat treatment NiTi wire technology.
- Greater resistance to cyclic fatigue.
- Improved file flexibility, while retaining cutting efficiency.
- SIE Sequence Kit includes a slider/shaper, and three finishing files.
- Ability to expand SIE sequence or create your own sequence to match the canal anatomy.

Reference: SIE

\$ 55<sup>00</sup>

Reference: A20


\$ 55<sup>00</sup>



**SCAN ME**  
For more information on  
F One

Product Name	21mm	25mm	31mm	Files per Pack
Style Italiano Endodontics Sequence #13/03 #20/04 #25/06 #35/04	SIE21	SIE25	-	3+1 Glide
#19/12 Orifice Opener 19mm	A20 2412 019 017	-	-	4 Open
#13/03 Shaper	SIE21GLIDE	SIE25GLIDE	-	4 Glide
#20/04 Small 4%	A20 1404 021 020	A20 1404 025 020	A20 1404 031 020	3+1 Open
#25/04 Primary 4%	A20 1404 021 025	A20 1404 025 025	A20 1404 031 025	3+1 Open
#30/04 Medium 4%	A20 1404 021 030	A20 1404 025 030	A20 1404 031 030	3+1 Open
#35/04 Large 4%	A20 1404 021 035	A20 1404 025 035	A20 1404 031 035	3+1 Open
#40/04 Extra Large 4%	A20 1404 021 040	A20 1404 025 040	A20 1404 031 040	3+1 Open
#20/06 Small 6%	A20 1406 021 020	A20 1406 025 020	A20 1406 031 020	3+1 Open
#25/06 Primary 6%	A20 1406 021 025	A20 1404 025 025	A20 1404 031 025	3+1 Open
#35/06 Large 6%	A20 1406 021 035	A20 1406 025 035	A20 1406 031 035	3+1 Open

## AF BLUE S4




Endo Motor settings:  
 Rotation: Fwd  
 Recommended Speed: 350rpm Max: 500rpm  
 Recommended Torque: 2 Nc Max Torque: 2.5Nc

Optional Files :  
 Glide Path: 13/02, 19/02  
 Working Length: 20/04, 25/04, 30/04, 40/04  
 20/06, 30/06, 35/06, 40/06



## AF BLUE Rotary



Endo Motor settings:  
 Rotation: Fwd  
 Recommended Speed: 350rpm Max: 500rpm  
 Recommended Torque: 2 Nc Max Torque: 2.5Nc

Optional Files :  
 Glide Path: 13/02, 19/02  
 Working Length: 20/04, 25/04, 30/04, 35/04, 40/04  
 20/06, 35/06



## AF Blue Rotary

*FLEXIBLE - EFFICIENT - FREEDOM OF CHOICE*  
*SIMPLE EASY TO FOLLOW SEQUENCES*

- Blue NiTi Wire Heat treatment technology.
- Greater resistance to cyclic fatigue.
- Improved file flexibility, while retaining cutting efficiency.
- S4 Sequence Kit includes an opener, slider/shaping, and two finishing files.
- Ability to expand S4 sequence or create your own sequence to match the canal anatomy.

Reference: S4 Sequence kits  
 (4 Files per pk)

**\$49<sup>50</sup>**

Reference: A04 (6 Files per pk)

**\$55<sup>00</sup>**



**SCAN ME**  
 For more information on  
 AF Blue

Product Name	21mm	25mm	31mm	Files per Pack
<b>S4 Sequence Kits</b> #17/12, #18/05, #25/06, #35/04	A04 1400 021 000	A04 1400 025 000	A04 1400 031 000	4
#19/12 Orifice Opener 19mm	A04 2612 019 017	-	-	6
#18/05 Shaper	A04 2605 021 018	A04 2605 025 018	-	6
#20/04 Small	A04 2604 021 020	A04 2604 025 020	-	6
#25/04 Primary	A04 2604 021 025	A04 2604 025 025	-	6
#30/04 Medium	A04 2604 021 030	A04 2604 025 030	-	6
#35/04 Large	A04 2604 021 035	A04 2604 025 035	-	6
#40/04 Extra Large	A04 2604 021 040	A04 2604 025 040	-	6
#20/06 Small 6%	A04 2606 021 020	A04 2606 025 020	-	6
#25/06 Primary 6%	A04 2606 021 025	A04 2606 025 025	-	6
#30/06 Large 6%	A04 2606 021 030	A04 2606 025 030	-	6
#35/06 Large 6%	A04 2606 021 035	A04 2606 025 035	-	6
#40/06 Extra Large 6%	A04 2606 021 040	A04 2606 025 040	-	6

# TX+™

## STERILE BLISTER PACKS

### Three Wire Technology in an easy to follow sequence.

- Silver Shape Memory Niti for TXA & TXPRG optimum control where needed.
- Gold Control Memory Niti for TX1 TX2 & TX3 for optimum flexibility to follow canal curvature and minimise canal transportation.
- Blue Control / Shape Niti for TX4 for a medium flexibility to encourage apical cutting efficiency and debris transportation.
- TX+ features with special geometry design creating a swaggering trajectory to enlarge debris removal efficiency.
- Variable taper design to have accurate shaping
- Existing workflow can be maintained




**SCAN ME**  
For more information on TX+

Reference: TX+


**Sterile**

**\$ 66<sup>00</sup>**

Product Name	21mm	25mm	31mm	Files per Pack
TXA, TXPRG, TX1, TX2, TX3, TX4 Assorted Pack	TXASSORT21	TXASSORT25	-	6
TXA Orifice Opener 19mm	TXA19	-	-	6
TXPRG Glider	TXPRG21	TXPRG25	-	6
TX1	TX1021	TX1025	TX1031	6
TX2	TX2021	TX2025	TX2031	6
TX3	TX3021	TX3025	TX3031	6
TX4	TX4021	TX4025	TX4031	6



### OPERATING STEPS




# VARIABLE TAPERED FILES

## VTaper GOLD

SIMPLE EASY TO FOLLOW SEQUENCE

Popular SX S1 S2 F1 F2 F3

- Gold Control Memory Niti for optimum flexibility to follow canal curvature and minimise transportation of the canal.
- Length marked
- Variable taper design for accurate shaping where needed
- Existing workflow can be maintained




**SCAN ME**  
For more information on  
VTaper Gold

Reference: V02

**\$55<sup>00</sup>**

Product Name	21mm	25mm	31mm	Files per Pack
SX, S1, S2, F1, F2, F3 Assorted Pack	V02 1600 021 000	V02 1600 025 000	V02 1600 031 000	6
SX Orifice Opener 19mm	V02 2604 019 019	-	-	6
S1	V02 2602 021 018	V02 2602 025 018	-	6
S2	V02 2604 021 020	V02 2604 025 020	-	6
F1	V02 2607 021 020	V02 2607 025 020	-	6
F2	V02 2608 021 025	V02 2608 025 025	-	6
F3	V02 2609 021 030	V02 2609 025 030	-	6




# W+™

STERILE BLISTER PACKS

FLEXIBLE, EFFICIENT, FAMILIAR SEQUENCE

- Radiopaque stopper
- Precise depth markings
- Reciprocating 30° 150° movement
- Variable regressive taper design to have accurate shaping
- Parallelogram cross section
- Safe guiding tip
- Existing workflow can be maintained




**SCAN ME**  
For more information on  
W+

Reference: W+


**Sterile**

**\$ 66<sup>00</sup>**

Product Name	21mm	25mm	31mm	Files per Pack
GLIDE, SML, PRIMARY x2, MED, LRG Assorted Pack	WASSORT21	WASSORT25	-	6
W+ Glider	WG021	WG025	-	6
W+ Small	WS021	WS025	WS031	6
W+ Primary	WP021	WP025	WP031	6
W+ Medium	WM021	WM025	WM031	6
W+ Large	WL021	WL025	WL031	6



## OPERATING STEPS



# RECIPROCATATING FILES

## AF BLUE R3

FLEXIBLE, EFFICIENT, FREEDOM OF SIZE CHOICE

- Reciprocating 30° 150° movement
- Blue NiTi Heat treatment technology for flexibility
- One File System
- Available in .04 & .06 Tapers
- Fixed Taper



**SCAN ME**  
For more information on  
AF Blue R3

Reference: A05

**\$ 55<sup>00</sup>**

Product Name	21mm	25mm	31mm	Files per Pack
#20/04 Small 4%	A05 2304 021 020	A05 2304 025 020	-	3
#25/04 Primary 4%	A05 2304 021 025	A05 2304 025 025	-	3
#20/06 Small 6%	A05 2306 021 020	A05 2306 025 020	A05 2306 031 020	3
#25/06 Primary 6%	A05 2306 021 025	A05 2306 025 025	A05 2306 031 025	3
#30/06 Medium 6%	A05 2306 021 030	A05 2306 025 030	A05 2306 031 030	3
#35/06 Large 6%	A05 2306 021 035	A05 2306 025 035	A05 2306 031 035	3
#40/06 Extra Large 6%	A05 2306 021 040	A05 2306 025 040	A05 2306 031 040	3

## R3 Tip Card

Irrigate, recapitulate with a #10 file and re irrigate after removing each rotary file

K File #10

Patency

13/02 16/02 19/02

Glide Path

350 rpm  
2N Torque

#20/06 SMALL

#25/06 PRIMARY 80% of cases

#30/06 MEDIUM

#35/06 LARGE

#40/06 EXTRA LARGE

Shaping and Finishing

150° Rev 30° Fwd  
300 - 350 rpm

- Scout the canal with hand use file
- Create glide path with your choice of 3 Glide path file sizes
- Choose proper size shaping files based on canal anatomy
- Shaping the canal with the chosen R3 File until it reaches working length
- After shaping ,finishing & Irrigation choose the matched GP points for obturation

## R+™

### RETREATMENT FILES


Reference: R+  
Size: 6 Files per Pack  
**Sterile** Blister Packs

**\$ 66<sup>00</sup>**



**SCAN ME**  
For more information on  
Retreatment Files


- Progressive on taper design and cutting flute length
- Remover of gutta-percha and other obturation material
- No. of Rings correspond to the three available sizes.
- Gold NiTi wire



D1 ISO 030, 9%, 16mm  
for coronal filling removal



D2 ISO 025, 8%, 18mm  
for middle filling removal



D3 ISO 020, 7%, 22mm  
for apical filling removal

Product Name	16mm	18mm	22mm	Files per Pack
R1	R116	-	-	6
R2	-	R218	-	6
R3	-	-	R322	6

## C PATH ROTARY


### SHAPE MEMORY GLIDE PATH WITH LASER LENGTH MARKS

Reference: F11  
Size: 6 Files per Pack

**\$ 55<sup>00</sup>**



**SCAN ME**  
For more information on  
Glide Path files



- Glide Path Management
- Respects the canal anatomy
- Rubber stopper marks the place of canal length and identify the direction of file
- Improved cyclic fatigue resistance
- Improved cutting efficiency
- Advanced inactive tip to avoid forming steps and to encourage debris upwards.
- Control memory Niti wire

Product Name	21mm	25mm	Files per Pack
#13/02	F21 2602 021 013	F21 2602 025 013	6
#16/02	F21 2602 021 016	F21 2602 025 016	6
#19/02	F21 2602 021 019	F21 2602 025 019	6

# SPECIALISED ROTARY FILES

## AF MAX

CLEAN - SHAPE - FINISH


Reference: F11  
Size: 6 Files per Pack

\$55<sup>00</sup>



**SCAN ME**  
For more information on  
AF Max Files

- Super elasticity and extreme flexibility
- Minimal stress applied to dentine walls
- Canal cleaning
- Excellent cleaning - finisher
- Remove smear layers efficiently
- GP removal possible
- Silver Niti wire



Product Name	21mm	25mm	Files per Pack
MAX 1 FINISHER #25/01	F11 2101 021 025	F11 2101 025 025	6
MAX 3 CLEAN & SHAPE #25/01	F11 2103 021 025	F11 2103 025 025	6

## AF CL


WORLDS ONLY LEDGE CORRECTION ROTARY FILE

Reference: F07  
Size: 6 Files per Pack

\$55<sup>00</sup>



**SCAN ME**  
For more information on  
Ledge Correction Files



- First ledge correction file in the world
- Pre curved tip
- The only specialist file available for ledged canals
- Silver NiTi wire

Product Name	21mm	25mm	Files per Pack
Assorted Pack			
#10 .06 - #15 .06 - #10 .08	F07 1600 021 000	F07 1600 025 000	6
#20 .07 - #25 .08 - #30 .09			
#10/06	F07 2606 021 010	F07 2606 025 010	6

## Hand Files

QUALITY K, C & H FILES



### K FILES

PERFECT BALANCE OF CUTTING EFFICIENCY  
& FLEXIBILITY

Square cross-section.

Stainless material. Length marked

Reference: KS  
Size: 6 Files per Pack  
**Sterile** Blister Packs

\$14<sup>00</sup>

### H FILES

HEDSTROEM FILES

Round cross-section for effective dentine cutting.

Stainless material. Length marked

Reference: HS  
Size: 6 Files per Pack  
**Sterile** Blister Packs

\$14<sup>00</sup>



### C FILES

STIFFER WIRE TECHNOLOGY

Specifically designed for narrow, calcified, and curved  
root canals. Heat-treated non cutting tip

Stainless material. Length marked

Reference: CS  
Size: 6 Files per Pack  
**Sterile** Blister Packs

\$18<sup>00</sup>



SCAN ME  
For more information on  
Hand Files



# STERILE HAND FILES

Product Name	21mm	25mm	31mm
Sterile K #06	KS 021 006	KS 025 006	KS 031 006
Sterile K #08	KS 021 008	KS 025 008	KS 031 008
Sterile K #10	KS 021 010	KS 025 010	KS 031 010
Sterile K #15	KS 021 015	KS 025 015	KS 031 015
Sterile K #20	KS 021 020	KS 025 020	KS 031 020
Sterile K #25	KS 021 025	KS 025 025	KS 031 025
Sterile K #30	KS 021 030	KS 025 030	KS 031 030
Sterile K #35	KS 021 035	KS 025 035	KS 031 035
Sterile K #40	KS 021 040	KS 025 040	KS 031 040
Sterile K #45	KS 021 045	KS 025 045	-
Sterile K #50	KS 021 050	KS 025 050	-
Sterile K #45-80	KS 021 458	KS 025 458	KS 031 458
Product Name	19mm	21mm	25mm
Sterile C #06	CS 019 006	CS 021 006	CS 025 006
Sterile C #08	CS 019 008	CS 021 008	CS 025 008
Sterile C #10	CS 019 010	CS 021 010	CS 025 010
Sterile C #12	CS 019 012	CS 021 012	CS 025 012
Sterile C #15	CS 019 015	CS 021 015	CS 025 015
Product Name	21mm	25mm	31mm
Sterile H #08	HS 021 008	HS 025 008	HS 031 008
Sterile H #10	HS 021 010	HS 025 010	HS 031 010
Sterile H #15	HS 021 015	HS 025 015	HS 031 015
Sterile H #20	HS 021 020	HS 025 020	HS 031 020
Sterile H #25	HS 021 025	HS 025 025	HS 031 025
Sterile H #30	HS 021 030	HS 025 030	HS 031 030
Sterile H #35	HS 021 035	HS 025 035	HS 031 035
Sterile H #40	HS 021 040	HS 025 040	HS 031 040
Sterile H #45	HS 021 045	HS 025 045	-
Sterile H #50	HS 021 050	HS 025 050	-
Sterile H #45-80	HS 021 458	HS 025 458	HS 031 458



**SCAN ME**  
For more information on  
Hand Files



## CERASEAL

### BIOCERAMIC SEALER

- Excellent sealing ability
- Moisture in the dentinal tubules and Calcium Silicate's chemical reaction
- Produces crystallization of Calcium Hydroxide
- CeraSeal is perfectly and fully hermetic in the root canal
- It prevents the propagation of bacteria
- CeraSeal never shrinks or expands in the root canal
- It prevents odontoclasia by keeping its stable volume
- The Single Cone technique obturation can be performed
- Shorter setting time
- Prevents from wash-out phenomenon

Product Name	Order Code	Price
CeraSeal	CERASEAL	\$198
10 Spare Tips	CERASEALTIPS	\$15

Reference: CERASEAL  
Size: 1x 2g Syringe + 15 Tips

\$198<sup>00</sup>



\$33<sup>00</sup>

## MD TEMP PLUS

### Hydraulic Temporary Material

- High resistance to drying out
- Easy application and removal
- Hermetic sealing of the cavity
- Excellent stability when set
- Proper crush strength
- White Colour in a 40g Jar



\$44<sup>00</sup>

## BINER LC

### Light Cured Cavity Liner

- Protect the pulp of the tooth
- Lining under all types of Fillings
- Protection when applying the total etch technique
- Light Curing type cavity liner, cures extremely hard
- Chemically bonds to adhesives and composites
- Insoluble in water and oral fluids
- Radiopacity



SCAN ME  
For more information on  
Sealers & Liners



## AD Seal Plus

RESIN BASED SEALER

- ADSEAL Plus contains Zirconium Oxide which does not cause tooth discolouration.
- ADSEAL Plus contains Calcium Phosphate, it offers biocompatibility and safety for humans.
- Provides the most appropriate flowability for RCT when base and catalyst are mixed.
- Shows about 25% higher radiopaque than competitor
- Offers a green-coloured dispense tip that enhances control of achieving the desired amount of sealer.
- Auto-mix tip option sold separately which provides users with easy and convenient mixing.

Product Name	Order Code	Price
AdSeal Plus	ADSEAL	\$77
20 Auto Mix Tips	ADSEALTIPS	\$22

Reference: ADSEALPLUS  
Size: 13.5g + Green Dispense Tip & Cap

Auto Mix Tips Sold Separate

\$77<sup>00</sup>



## Opacal

CALCIUM HYDROXIDE DRESSING

Opacal is a radiopaque, creamy and ready-to-use calcium hydroxide paste with a pH above 12  
Paste delivery in root canal is ensured by dedicated application tip

### Composition

Calcium Hydroxide 30% Barium Sulfate 24%

### Key benefits

- High pH value
- High radiopacity
- Creamy consistency
- Easy and comfortable paste delivery
- Easy to remove (non-setting material)
- Rapid ionic release

Product Name	Order Code	Price
OPACAL	OPACAL	\$77
10 Spare Tips	CERASEALTIPS	\$15

Reference: OPACAL  
Size: 1x 3ml Syringe + 10 Tips

\$77



SCAN ME  
For more information on  
Sealers & Dressings



### High Precision

#### Gutta Percha Points

Length Marked  
Matched Tapers & Sizes

Reference: MMGP  
Size: 60 Per Pack

\$ 19<sup>80</sup>

4% Taper	6% Taper	V-Taper
#20	#20	F1
#25	#25	F2
#30	#30	F3
#35	#35	-
#40	#40	F1 - F2 - F3



### High Absorption Paper Points

Length Marked  
Matched Tapers & Sizes

Reference: PP  
Size: 60 Per Pack

\$ 17<sup>60</sup>

4% Taper	6% Taper	V-Taper
#20	#20	F1
#25	#25	F2
#30	#30	F3
#35	#35	-
#40	#40	F1 - F2 - F3

# GUTTA PERCHA & PAPER POINTS



## Lentulos - Paste Carrier


Spiral Paste Fillers are made from high quality stainless steel. They are also autoclavable and colour coded. Used for distribution of pastes and medicaments into the prepared canals.

Sizes #25, #30, #35, #40

Length:25mm



Size	Order Code	Description
#25	PF2525	25mm Pack of 4
#30	PF2530	25mm Pack of 4
#35	PF2535	25mm Pack of 4
#40	PF2540	25mm Pack of 4
#25-#40	PF25A	25mm Pack of 4



## Gates Glidden

High quality medical Stainless Steel used for expanding the canal entrance or preparation for middle-upper part of the canals.

Short blades, the diameter for working blades is from 0.50mm - 1.50mm.

Ring numbers show different sizes.


Recommended speed: 800-1200 rpm.

Sizes:#1, #2, #3, #4, #5, #6

Length: 28mm



Size	Order Code	Description
#1	GTG2801	28mm Pack of 6
#2	GTG2802	28mm Pack of 6
#3	GTG2803	28mm Pack of 6
#4	GTG2804	28mm Pack of 6
#5	GTG2805	28mm Pack of 6
#6	GTG2806	28mm Pack of 6



## Peeso Reamers

High Quality Medical Stainless Steel for preparing the 1/3 coronal part of canals, obturation, post removal and in preparation for posts.

Long blades, diameter for working blades is from 0.70mm-1.70mm.

Ring numbers show different sizes.

Recommended speed: 800-1200 rpm.

Sizes: #1, #2, #3, #4, #5, #6

Length:28mm



Size	Order Code	Description
#1	PSR2801	28mm Pack of 6
#2	PSR2802	28mm Pack of 6
#3	PSR2803	28mm Pack of 6
#4	PSR2804	28mm Pack of 6
#5	PSR2805	28mm Pack of 6
#6	PSR2806	28mm Pack of 6



## RetroMTA®

Used for vital pulp therapy.  
Pulp capping, pulpotomy and pulpectomy in vivo condition.

Initial Ph is up to 12.5 and down to 7.8~8.0 in 4 weeks

- Vital Pulp Therapy & Root Repair
- Fast setting
- No discoloration
- No heavy metal
- No cell toxicity
- Good handling properties

Reference: 3006 40 2000

Size: 0.3g x 2 Caps

Reference: 3006 40 2001

Size: 0.3g x 8 Caps

\$ 55<sup>00</sup>

\$ 176<sup>00</sup>



RetroMTA is a hydraulic Bioceramic material for repair and vital pulp therapy.

It has no discoloration even in case of blood contamination.



## OrthoMTA®

- Easy and fast root repair
- Initial setting of 3 hours and final setting of 6 hours
- No heavy metals
- No expansion(0.09%)
- Removable
- Economic package and price

Reference: 3006 40 2002

Size: 0.2g x 5 Vials

Reference: 3006 40 2003

Size: 0.2g x 20 Vials

\$ 99<sup>00</sup>


\$ 330<sup>00</sup>



OrthoMTA is an innovative orthograde root canal grafting material that is unaffected by moisture or blood contamination. Ideal synthetic MTA for orthograde root canal filling. High purity MTA without hexavalent chromium, a carcinogenic heavy metal, and does not shrink after root canal grafting / filling.



**SCAN ME**  
For more information on  
MTA Material



## Bio MTA

### MTA CARRIER SYSTEM

Provides unrestricted access to the canal orifice, allowing an accurate orthograde grafting procedure with MTA whilst avoiding the dispersion of the root grafting material in the oral cavity.

The OrthoMTA grafting system means

1. Reduced manual work
2. No micro-leakage
3. No waste of material
4. Increased positive prognosis

#### Items included in the Biofilling Starter Kit

- **OrthoMTA Carrier**
- **4 vials of OrthoMTA**
- **Plunger OMP-03**
- **Bender**
- **Plunger 18g, 16g**
- **Needle Tip 18g, Needle Tip 16g**
- **Resin Block**

Product Name	Reference	Price
Biofilling Starter Kit	9018 49 8000	\$990
Disposable needle 16g x 50	9018 49 8005	\$88
Disposable Needle 18g x 50	9018 49 8003	\$88
Plunger 16g x 4	9018 49 8004	\$165
Plunger 18g x 4	9018 49 8002	\$165
Needle Bender x 3	9018 49 8006	\$22
Replacement Carrier	9018 49 8001	\$385

Reference: 9018 49 8000

\$ 990<sup>00</sup>



**SCAN ME**  
For more information on  
MTA Carrier



## NEXCORE

DUAL CURED CORE BUILD UP MATERIAL

2 Shades - Tooth or White for ceramic restoration

Post Adhesions and core build-up in restorations

Reference: NEXCOREKIT

Colour: White

Colour: Tooth

\$99<sup>00</sup>



- Suitable viscosity and flowability for post cementation and core build-up
- Dual Curing: Light Curing and self curing are possible, stress free from polymerisation failure of resin
- Biocompatibility (ISO10993)
- Easy to mix by auto mixing, direct injection intraoral and endo type tips
- Radiopaque
- 2 Shades - Tooth, White for ceramic restoration
- Similar cuttability like dentine (Flexural modulus : NexCore 14.2GPa)



### DUAL SYRINGE

2 x 9g Syringes for Post adhesions and Core Build-up in restorations



### MIXING TIPS

20 x Easy auto mixing tips



### ENDO TIPS - INTRORAL TIPS

10 x Direct injection Intraoral  
10 x Endo type tips



**SCAN ME**  
For more information on  
Core Build up



# SYRINGES, NEEDLES & DAM



## Syringes

Various Endodontic Irrigation needles

3ml, 5ml, 10ml Luer Lock or Curved tip utility syringes

Product Name	Reference	Price
3ml Luer Sterile x 100	Z0037	\$24.20
5ml Luer Sterile x 100	Z0038	\$30.25
10ml Luer Sterile x 100	Z0039	\$36.90
Utility x 50	49001	\$44



## Irrigation Needles

Safe End Dual Port Irrigation

30G Closed End, Dual side vent.  
Even distribution of irrigants, reduced apical pressure. Sterilised individual packs with capped tips.

Reference: 1120024  
Size: 30g x 25mm

\$99<sup>00</sup>



## Latex Free Rubber Dam

Give a Dam !

Powder Free, Latex Free.


Medium - Mint

Silky smooth rubber dam in a soft blue giving excellent eye relief.

\$44<sup>00</sup>

Reference: GIVEADAM

Size: 15 pieces - 152mm x 152mm



## Tapered Irrigation Needles

Safe End Dual Port Irrigation

Unlike conventional metal needles, IrriFlex® easily adapts to complex root canal anatomies. The flexibility of the polypropylene body allows the needle to access the apical region without resistance or damage to the dentinal walls.

Reference: PDI20  
Size: 20 x Individually wrapped

\$74<sup>80</sup>



SCAN ME  
For more information on  
Irrigation

\$ 2350<sup>00</sup>



**E Connect S +**  
In Built Apex Locator

Cordless, Brushless motor.  
All Files systems including Forward and Reciprocation.  
Fully adjustable in all settings.  
OLED Colour LCD Display.  
Torque Calibration.  
Ambidextrous Display.

Product Name	Price
E Connect S+ Complete Kit including 1x Contra Angle Head	\$ 2350
Spare Parts	Price
Contra Angle 1500rpm Head	\$550
White USB C Apex Cable	\$44
Black Apex File Holder	\$44
Lip Clip	\$16.5

Different working modes



- CW: Clockwise Continuous Rotation
- CCW Counter Clockwise Continuous Rotation
- REC: Reciprocation Movement
- ATC: Adaptive Torque Control
- Auto Start and Stop
- Apical Reverse
- Apical Slow Down



Technical

- Rotational speed range: 120 - 1500 rpm
- Torque range: 0.5 - 5.0 Ncm
- Memory programmes: 10
- Powerful Li-ion battery 1900 mAh



**SCAN ME**  
For more information on  
Endo Motors



**\$ 1150<sup>00</sup>**



## Extreme

Slim Light Egonomic

Cordless, Brushless motor.

All Files systems including Forward and 5 preset Reciprocation settings.

OLED Colour LCD Display.

Torque Calibration.

Ambidextrous Display.

Product Name	Price
Extreme Complete Kit including 1x Contra Angle Head	\$ 1150
Spare Parts	Price
Contra Angle 650rpm Head	\$330
White USB C Apex Cable	\$44
Black Apex File Holder	\$44
Lip Clip	\$16.5



### Different working modes

- CW: Clockwise Continuous Rotation
- CCW Counter Clockwise Continuous Rotation
- REC: Reciprocation Movement



### Technical

- Rotational speed range: 120 - 650 rpm
- Torque range: 0.5 - 4.0 Ncm
- Memory programmes: 10
- Powerful Li-ion battery 800 mAh



**SCAN ME**  
For more information on  
Endo Motors

\$ 1995<sup>00</sup>



## E Connect S

In Built Apex Locator



- Cordless, Brushless motor.
- All Files systems including Forward and Reciprocation.
- Fully adjustable in all settings.
- OLED LCD Display.
- Torque Calibration.
- Ambidextrous Display.

Product Name	Price
E Connect S+ Complete Kit including 1x Contra Angle Head	\$ 1995
Spare Parts	Price
Contra Angle Head	\$550
White Micro USB Apex Cable	\$44
Black Apex File Holder	\$44
Lip Clip	\$16.5

### Different working modes



- CW: Clockwise Continuous Rotation
- CCW Counter Clockwise Continuous Rotation
- REC: Reciprocation Movement
- ATC: Adaptive Torque Control
- Auto Start and Stop
- Apical Reverse
- Apical Slow Down



### Technical

- Rotational speed range: 120 - 1000 rpm
- Torque range: 0.5 - 4.0 Ncm
- Memory programmes: 10
- Powerful Li-ion battery 1500 mAh



**SCAN ME**  
For more information on  
Endo Motors



\$880.00

## Ultra X

Irrigant Activator



- Disrupt the smear layer and bio-film and open up the plugged dentinal tubules
- Tip oscillation and vibration in frequency of 45KHz to amplify the efficiency of irrigants
- GP removal during retreatment
- Broken instrument retrieval

Product Name	Price
Ultra X Complete Kit including 2 of each Colour Tips (6 Tips total)	\$ 880
Spare Parts	Price
Silver Tips 2x 20/02 21mm	\$55
Blue Tips 1x 20/02 18mm 1x 25/02 18mm	\$55
Gold Tips 1x 20/02 18mm 1x 25/02 18mm	\$55

### Irrigation Protocol Using Ultra X



- Prepare the root canal to the desired shape
- Use a side-vented needle to inject the final rinse solution into the canal
- Cover the device with a protective sleeve to avoid surface damages



- Select an activator tip that fits passively when placed 2-3 mm short of working length

- Agitate the solution using short vertical strokes for around 30 seconds.
- Be careful to not breach the apical 2 mm
- Replenish the irrigant and use suction to remove loose debris
- Repeat the agitation and replenishment cycle 4-5 times
- The wrench, protective sleeves, and tips are re usable and autoclavable.




**SCAN ME**  
For more information on  
Ultra X



\$ 440<sup>00</sup>

## E Pex Pro

Smart, Accurate and Reliable Apex Locator



- Advanced Multi-frequency Technology
- High Precision in Wet & Dry Canals
- Large 3.5" LCD Display
- Automatic Calibration
- Compact Size & Stable Design
- Powerful 1500mAh Lithium-ion Battery

Product Name	Price
EPex Pro Complete Kit including 1x Micro USB Cable 2x Lip Clips & 2x File Holders	\$ 440
Spare Parts	Price
White Micro USB Apex Cable	\$44
Black Apex File Holder	\$44
Lip Clip	\$16.5

### FEATURES



- 4th Generation - Highly developed multi-frequency based network impedance measurement technology.
- LCD Display -Bright and clear display.
- Apical foramen of High precision localization in wet and dry canals.



- Automatic calibration - Ensures steady and correct measurements.
- Stable design and compact size allows for outstanding suppleness during procedures.



**SCAN ME**  
For more information on  
Apex Locators



\$ 3960<sup>00</sup>



## Hyperlight

Low patient dosage - TGA Approved

Smart Angulation system  
 Lightweight, cordless design  
 One - handed operation possible  
 Long lasting battery. 500 X-Ray shots on a single battery charge  
 Excellent image quality with 0.4mm focal spot  
 Internal leaded double shield protects the operator

Low patient dosage, TGA & FDA approved

Our Hyperlight X-Ray machine is compact in size and its portable nature helps you to use it in a comfortable manner. Moreover, this unit is highly ensured with safety standards owing to its minimum exposure time.

Fully adjustable setting to be used with all Xray sensors and Plate systems..

Product Name	Price
Hyperlight Complete Kit including 1x Protection Shield	\$ 3960

### Specifications



- Model: Hyper light
- Focal Spot: 0.4mm
- Tube Voltage: 65kV
- Tube Current: 2.5mA
- Weight: 1.8 kg
- Battery: : 2500mAh x 4
- Exposure Time Setting: 0.02-2 Sec



### Packaging

- 1 x HyperLight Portable DC X-Ray Machine
- 1x Charger
- 1 x Protection shield



**SCAN ME**  
 For more information on  
 Hyperlight

\$ 462<sup>00</sup>



### Curing Pen

Sleek Light Cure Device

Three Curing Option

1. Ramp, Pulse
2. Caries
3. Detection Mode

Caries detection mode: 600mw/cm<sup>2</sup>

Working light intensity 1000 mw/cm<sup>2</sup>, 1500 mw/cm<sup>2</sup>, 2300 mw/cm<sup>2</sup>

4 LEDs curing light - 3 Blue & 1 Violet

Wavelength range: 385 - 515 nm

360 rotatable head

Wide Digital display

Consistent light output

Product Name	Price
CurePen Complete Kit including 1x Head and 100x Sleeves	\$ 462
100x Barrier Sleeves	\$ 24.2



#### Specifications

- Light Weight: 109 gm
- Ergonomic design
- Dimensions 24\*21\*230 mm
- Battery: >1200 times (3s at max intensity)




#### Packaging

- 1 x Curing Pen with Head
- 1x USB C Charger + Base
- 100 x Barrier Sleeves



**SCAN ME**  
For more information on  
Curing Lights





**\$1320<sup>00</sup>**

\*Satelec and EMS are trademarks of Acteon® and EMS®, respectively.

## Ultra Mint Pro

Scaling | Periodontic | Endodontic

- Wide range of frequency: 25kHz-42kHz
- EMS and Scaler type\* tip threads are compatible with UltraMint & UltraMint Pro
- 900ml liquid container, can work directly with sodium hypochlorite, chlorhexidine, hydrogen peroxide and saline
- Internal self-cleaning container.
- Detachable and autoclavable handpiece holder, resist cross- infection
- Ergonomic foot pedal
- Water control directly from the handpiece
- Intelligent negative Feedback, In accordance to different tooth surfaces, it will automatically adjust the pattern to the tip, offering a stable and sustained scaling treatment that is comfortable and painless for the patient.

Product Name	Price
Scaler System with Water Tank including 1x Handpiece, 6x Tips, 1x Wrench	\$ 1320
Scaler System without Water Tank including 1x Handpiece, 6x Tips, 1x Wrench	\$ 990
Spare Parts	
Satalec* type Hand Piece Kit with tips & Wrench	\$440
EMS* type Hand Piece Kit with tips & Wrench	\$440
Satalec* type Hand Piece only	\$231
EMS* type Hand Piece only	\$231
Ultra Mint Torque Wrench	\$88
Ultra Mint Water Tank & Cover	\$88

### Inclusions



- Scaler unit with 900ml container
- Detachable handpiece LED-autoclavable
- Detachable and autoclavable handpiece holder
- Foot pedal
- Autoclavable plastic box
- Comes with 6 x Tips



- Choose either, (EMS type\* tip thread is sent as default, if you want Satelec\* type tip thread please advise.
- EMS type tips included: G1,G2,G4,P1,P3,E4.
- Satelec type tips included: E3-SD, E5-S, E7-S, P4-S, G1-S, G2-S
- Universal wrench

\*Satelec & EMS are trademarks of their respective owners



**SCAN ME**  
For more information on  
Ultra Mint Pro



Fast Fill 23g Tip  
 Size: 23g  
 Qty: 1 per pack

**\$ 77<sup>00</sup>**



Fast Fill 25g Tip  
 Size: 25g  
 Qty: 1 per pack

**\$ 77<sup>00</sup>**



GP Pluggers  
 Size: #35/70, #40/80, #50/100, #60/120  
 Qty: Set of Four

**\$ 175<sup>00</sup>**



Fast Fill GP Bars  
 Size: -  
 Qty: 100 per pack

**\$ 66<sup>00</sup>**



**Star Shape  
 Fast Pack Pro**

Fast Pack Pro Tips  
 Size: Black or Yellow or Blue  
 Qty: 1 per pack

**\$ 77<sup>00</sup>**



**Hex Shape  
 Fast Pack**

Fast Pack Tips  
 Size: Black or Yellow or Blue  
 Qty: 1 per pack

**\$ 77<sup>00</sup>**

\$ 1100<sup>00</sup>



\$ 1540<sup>00</sup>



## Fast Fill

Back Fill for 3D Obturation System

- Motorized hand-piece allows for a consistent flow and homogenous 3D filling
- Pencil-grip hand-piece enhances control
- Precise temperature control for consistency
- Extra-long needles offer convenience and unprecedented visibility
- Hand-piece shows remaining GP
- Heats quickly

Product Name	Price
Fast Fill Complete Kit including 2x Tips 100 GP Pellets	\$ 1540



## Fast Pack Pro

Down pack for 3D Obturation System

- Only the working end of the tip is heated- avoids need for insulated sleeve
- Non-working end stays cool- making it comfortable for patients
- 360 degree start/stop ring, allows operation from any angle comfortable
- 4 hours of continuous operation at 250 degree Celsius
- 2 hours to complete charge
- 15 minutes turbo-charge to 30%
- 50 cases with a single charge

Product Name	Price
Fast Pack Pro Complete Kit including 3x Heating Tips	\$ 1100



**SCAN ME**  
For more information on  
Obturation



# CMS One Step Obturation System

## HANDLELESS CARRIER

Reference: OSDTR2060  
1x Oven, 2x Tweezers,  
12x Verifiers, 20x Assorted GP

\$440<sup>00</sup>



OneStep® is the second generation of obturators from CMS Dental. The product was developed due to a request from dentists who wanted an obturator without a handle. It gives a better working area and space to perform the treatment.



OneStep® obturator consists of a radiopaque plastic core surrounded by natural thermo-plastic gutta-percha. The plastic core is available in ISO sizes # 20-60.



### FAST

Warming up an obturator takes 1 min.  
Up to 4 obturators can be hung in the oven at a time.



### EASY

Hold the tweezers so that the tip of the tweezers corresponds to the working length before placing it in the oven.  
Then simply insert into the prepared root canal.



### SAFE

More than 5 million root fillings have been made with One-Step which has proven to be a quality product.



## OneStep®

A SAFE AND EFFECTIVE ROOT FILLING SYSTEM


OneStep® products are made with natural gutta-percha, a material that has been used for root filling for well over 100 years.

Root canal filling with OneStep® is easy and provides a perfect apical seal.


OneStep® can be used for root canal fillings instrumented with traditional hand files as well as with all known rotary file systems on the market.

- **One Step is an easy and safe treatment method that provides a dense root filling**
- **The solid core ensures:**
  - a) **That gutta-percha gets right down to the apex.**
  - b) **That a certain lateral pressure is applied, forcing the gutta-percha into lateral canals and outlets.**
- **A control X-ray should always be taken to confirm a correct and dense root filling.**

Product Name	Reference	Price
OneStep Complete Kit 1x Oven, 20x GP Obturators, 2x Tweezers & 12x Verifiers	OSDTR2060	\$ 440
Consumables	Reference	Price
GP #20 Refill Pack of 20	OSR2020	\$ 124
GP #25 Refill Pack of 20	OSR2025	\$ 124
GP #30 Refill Pack of 20	OSR2030	\$ 124
GP #35 Refill Pack of 20	OSR2035	\$ 124
GP #40 Refill Pack of 20	OSR2040	\$ 124
GP #50 Refill Pack of 20	OSR2050	\$ 124
GP #60 Refill Pack of 20	OSR2060	\$ 124
Verifier #20 Pack of 6	SSV620	\$ 24.9
Verifier #25 Pack of 6	SSV625	\$ 24.9
Verifier #30 Pack of 6	SSV630	\$ 24.9
Verifier #35 Pack of 6	SSV635	\$ 24.9
Verifier #40 Pack of 6	SSV640	\$ 24.9
Verifier #50 Pack of 6	SSV650	\$ 24.9
Verifier #60 Pack of 6	SSV660	\$ 24.9
Locking Tweezer 100°	OSOTWE100N	\$ 88
Locking Tweezer 135°	OSOTWE135N	\$ 88



**SCAN ME**  
For more information on  
Obturation



### Munce Discovery Burs® Troughers

The Shallow Troughers (31mm) were specifically designed for procedures on the chamber floor or just beyond such as troughing for the MB2 canal of maxillary molars, troughing for the middle mesial canal of mandibular molars or to “brush out” the isthmus between the MB & ML canals of mandibular molars. Shallow Troughers are also useful where limited inter-arch space precludes the use of the Deep Troughers.

The Deep Troughers (34mm) were designed for deeper operations within the root such as exposure of separated instruments, removal of pulpal tissue from deep within C-shaped canals and chasing severely-calcified canals into the mid-root and beyond.

The Cariesectomy Bur (31mm) has a #6 round carbide head and was designed to facilitate the type of deep caries removal often required in endodontic cases, and for which there was previously no comparable instrument. This bur—often in combination with the #4 Shallow Trougher—is ideal for caries removal from around crown margins from within the access cavity.

- Features
- Slow-speed
  - Stiff Shaft
  - No Heat Generation
  - Maintains a Functional View Corridor beyond the Handpiece Head
  - Latch-type
  - Durable Carbide Head

Product Name	31mm	34mm	Price
Variety pack 6 Pieces/Pk	S0VP	D0VP	\$ 137.5
# ¼ Pink 4 Pieces/Pk	S4025	D4025	\$ 99
# ½ Grey 4 Pieces/Pk	S4050	D4050	\$ 99
# 1 Purple 4 Pieces/Pk	S4100	D4100	\$ 99
# 2 White 4 Pieces/Pk	S4200	D4200	\$ 99
# 3 Yellow 4 Pieces/Pk	S4300	D4300	\$ 99
# 4 Red 4 Pieces/Pk	S4400	D4400	\$ 99
# 6 Green 4 Pieces/Pk	S4600	-	\$ 99

**\$137<sup>50</sup>**

## Comprehensive Endo Bur Kit



Shape	Code	Size	Application
	A4PF		Safe ceramic penetration without microfractures
	FGB331		Fast & easy metal penetration
	ENDOZ		Pulp chamber opening without risk
	E3RM		Metal post removal / meta screw weakening cementation
	E3RF		Safe search of canal entrance / span
	A1LM		Finding the canal entrance / span
	A4LM		Pulp chamber opening
	PR2SF		Middling composite polishing
	PR2XF		Gentle composite polishing & finishing
	A4SF		Middling composite polishing
	A4XF		Gentle composite polishing & finishing
			11 Pieces + Aluminium Bur Block

**Buy 10  
Get 10%  
OFF**  
Mix n Match all  
Burs



**\$36<sup>30</sup>**

### Endo Z Safe End

Code	Size	Price
ENDOZ	016 Non cutting Tip	\$ 36.3
3 Pieces per Pack		



**\$24<sup>20</sup>**

### B1558 Metal cutting

Code	Size	Price
B1558	012 Carbide	\$ 24.2
3 Pieces per Pack		



**\$24<sup>20</sup>**

### FGB Metal Cutting

Code	Size	Price
FGB330	008 Carbide	\$ 24.2
FGB331	010 Carbide	\$ 24.2
3 Pieces per Pack		



**\$24<sup>20</sup>**

### FG4OS Carbide Rose

Code	Size	Price
FG4OS	014 Carbide	\$ 24.2
3 Pieces per Pack		



**SCAN ME**  
For more information on  
Burs



\$ 24<sup>20</sup>

### Extra Long Surgical

Code	Size	Price
A2XLM	Medium 014	\$ 24.2
A3XLM	Medium 018	\$ 24.2
2 Pieces per Pack		



\$ 33<sup>00</sup>

### Round

Code	Size	Price
A4PF	Magic Touch 022 2/Pk	\$ 24.2
A4SF	Super Fine 022 6/Pk	\$ 24.2
A4XXF	Extra Extra Fine 022 3/Pk	\$ 33



\$ 19<sup>80</sup>

### Separation

Code	Size	Price
E3RF	Fine 008	\$ 19.8
E3RM	Medium 008	\$ 19.8
6 Pieces per Pack		



\$ 19<sup>80</sup>

### Round End Taper

Code	Size	Price
F1RXLF	Fine 016	\$ 19.8
F1RXLM	Medium 016	\$ 19.8
2 Pieces per Pack		



\$ 33<sup>00</sup>

### Pear

Code	Size	Price
PR2SF	Super Fine 018 6/Pk	\$ 33
PR2XF	Extra Extra Fine 018 3/Pk	\$ 33



\$ 24<sup>20</sup>

### Endo Access

Code	Size	Price
PR1210C	Coarse 012/016	\$ 24.2
PR1910C	Coarse 018/019	\$ 24.2
2 Pieces per Pack		



\$ 24<sup>20</sup>

### Endo Z Safe end Diamond

Code	Size	Price
L4M	Medium 016	\$ 24.2
2 Pieces per Pack		



\$ 19<sup>80</sup>

### Round end Taper

Code	Size	Price
C0C	Coarse 012	\$ 19.8
C1C	Coarse 014	\$ 19.8
M3M	Medium 018	\$ 19.8
6 Pieces per Pack		



\$ 24<sup>20</sup>

### Round Minimally Invasive

Code	Size	Price
MA0M	Medium 009	\$ 24.2
MA1M	Medium 012	\$ 24.2
MA2M	Medium 014	\$ 24.2
MA25M	Medium 016	\$ 24.2
3 Pieces per Pack		



\$ 24<sup>20</sup>

### Tapered Minimally Invasive

Code	Size	Price
MAB0M	Medium 008	\$ 24.2
MAB3M	Medium 010	\$ 24.2
3 Pieces per Pack		



\$ 19<sup>80</sup>

### Extra Long Round end Taper

Code	Size	Price
F1RXLF	Fine 016	\$ 19.8
F1RXLM	Medium 016	\$ 19.8
2 Pieces per Pack		



\$ 44<sup>00</sup>

### Surgical Length Round

Code	Size	Price
A1LM	Medium 012	\$ 44
A2LM	Medium 014	\$ 44
A2LC	Coarse 014	\$ 44
A3LM	Medium 018	\$ 44
A3LC	Coarse 018	\$ 44
A4LM	Medium 022	\$ 44
6 Pieces per Pack		



# Diamond Burs & Instruments

## TWEEZERS



27121

\$44<sup>00</sup>

Reference	Description
27121	Tweezer College Self Locking

## PROBES



27550

\$44<sup>00</sup>



Round: Solid Handle - Double Ended

Reference	Description
27550	Probe Endodontic DG 16 D/E

## EXCAVATORS



27540N

\$52<sup>80</sup>



27565N

\$66<sup>00</sup>



'Lite' 10mm Hollow Handle

Reference	Description
27540N	Excavator D/E 31L
27565N	Glick Endospoon #2

## SPREADERS



27551



27552



27553



27560



27561

\$44<sup>00</sup>



Round: Solid Handle - Single Ended

Reference	Description
27551	Spreader S/E 40
27552	Spreader S/E 60
27553	Spreader S/E 80
27560	Spreader D/11
27561	Spreader D/11T

## MIRRORS



27100

\$15<sup>40</sup>



60416

\$44<sup>00</sup>

Reference	Description
27100	Mirror Simple Stem Handle Round
60416	Mirror Examination Front Surface Plane #4




SCAN ME  
For more information on  
Instruments

EMS G1  
Satelec G1-S

\$1650

Function Description: All-round removal of supragingival calculus



Mode: G  
Power: 1-10


Recommended Service Life  
200 Teeth

Water Flow Rate  
mid&high

EMS G2  
Satelec G2-S

\$1650

Function Description: Removal of significant supragingival calculus



Mode: G  
Power: 1-10


Recommended Service Life  
200 Teeth

Water Flow Rate  
mid&high

EMS G4  
Satelec G4-S

\$1650

Function Description: Used for supragingival and interdental space scaling



Mode: G  
Power: 1-10


Recommended Service Life  
200 Teeth

Water Flow Rate  
mid&high

EMS G5  
Satelec G5-S

\$1650

Function Description: Removal of significant supragingival calculus



Mode: G  
Power: 1-10


Recommended Service Life  
200 Teeth

Water Flow Rate  
mid&high

EMS P1  
Satelec P1-S

\$2200

Function Description: Removal of subgingival calculus and periodontal irrigation



Mode: P  
Power: 1-10


Recommended Service Life  
150 Teeth

Water Flow Rate  
mid&high

EMS P3  
Satelec P3-S

\$2200

Function Description: Removal of calculus deep in periodontal pocket and its irrigation



Mode: P  
Power: 1-6


Recommended Service Life  
150 Teeth

Water Flow Rate  
mid&high

EMS P4  
Satelec P4-S

\$2200

Function Description: Removal of calculus in shallow pockets



Mode: P  
Power: 1-6


Recommended Service Life  
200 Teeth

Water Flow Rate  
mid&high

EMS P52  
Satelec P52-S

\$3300

Function Description: Supragingival scaling



Mode: P  
Power: 1-10


Recommended Service Life  
150 Teeth

Water Flow Rate  
high

EMS P50L  
Satelec P50L-S

\$3300

Function Description: Used for subgingival calculus probing, subgingival proximal surface scaling and root planing



Mode: P  
Power: 1-10


Recommended Service Life  
150 Teeth

Water Flow Rate  
high

EMS P50R  
Satelec P50R-S

\$3300

Function Description: Used for subgingival calculus probing, subgingival proximal surface scaling and root planing



Mode: P  
Power: 1-10


Recommended Service Life  
150 Teeth

Water Flow Rate  
high

EMS P53L  
Satelec P53L-S

\$3300

Function Description: Used for subgingival calculus probing, subgingival proximal surface scaling and root planing



Mode: P  
Power: 1-10


Recommended Service Life  
150 Teeth

Water Flow Rate  
mid&high

EMS P53R  
Satelec P53R-S

\$3300

Function Description: Used for subgingival calculus probing, subgingival proximal surface scaling and root planing



Mode: P  
Power: 1-10


Recommended Service Life  
150 Teeth

Water Flow Rate  
mid&high

EMS P59  
Satelec P59-S

\$3300

Function Description: Used for subgingival calculus probing and scaling



Mode: P  
Power: 1-6

Recommended Service Life  
150 Teeth


Water Flow Rate  
high

EMS P3D  
Satelec P3-SD

Diamond Coated

\$2750

Function Description: Used to level off the surface of endo after the periodontal flap surgery



Mode: P  
Power: 1-6

Recommended Service Life  
50 Teeth


Water Flow Rate  
mid&high

EMS P4D  
Satelec P4-SD

Diamond Coated

\$2750

Function Description: For locating root canal and removal of calcification on the one-third of crown



Mode: P  
Power: 1-6


Recommended Service Life  
50 Teeth

Water Flow Rate  
low, mid&high

EMS E3  
Satelec E3-S

\$2200

Function Description: For lateral condensation of gutta percha by heating effect



Mode: E  
Power: 1-3

Recommended Service Life  
500 root canals

Water Flow Rate  
No water




SCAN ME  
For more information on  
Scaler Tips

EMS E4  
Satelec E4-S

\$2750

Function Description: Removal of foreign filling material during re-treatment



Mode: E  
Power: 1-6


Recommended Service Life  
100 root canals

Water Flow Rate  
No water

EMS E5  
Satelec E5-S

\$2750

Function Description: Removal of foreign filling material during re-treatment



Mode: E  
Power: 1-6


Recommended Service Life  
100 root canals

Water Flow Rate  
No water

EMS E7  
Satelec E7-S

\$2200

Function Description: Removal of broken instrument



Mode: E  
Power: 1-3


Recommended Service Life  
50 root canals

Water Flow Rate  
mid&high

EMS E9  
Satelec E9-S

\$2750

Function Description: Flexible, used for irrigation



Mode: E  
Power: 1-3

Recommended Service Life  
500 root canals


Water Flow Rate  
No water

EMS E3D  
Satelec E3-SD

Diamond Coated

\$3300

Function Description: Enlargement of the canal walls to remove the broken instrument



Mode: E  
Power: 1-3


Recommended Service Life  
50 root canals

Water Flow Rate  
low,mid&high

EMS E30  
Satelec E30-S

\$2750

Function Description: Removal of deep foreign particles in canal



Mode: E  
Power: 1-3


Recommended Service Life  
50 root canals

Water Flow Rate  
low&mid

EMS E31  
Satelec E31-S

\$2750

Function Description: Removal of foreign particles in one-third of the canal



Mode: E  
Power: 1-3

Recommended Service Life  
50 root canals


Water Flow Rate  
low&mid

EMS E4D  
Satelec E4-SD

Diamond Coated

\$3300

Function Description: Removal of hard material in canal during re-treatment



Mode: E  
Power: 1-3

Recommended Service Life  
50 root canals

Water Flow Rate  
No water

**EMS E5D**  
**Satelec E5-SD**

**Diamond Coated**

**\$3300**

Function Description: Removal of hard material in canal during re-treatment

Mode: E  
Power: 1-3

Recommended Service Life  
50 root canals

Water Flow Rate  
No water

**EMS E6D**  
**Satelec E6-SD**

**Diamond Coated**

**\$3300**

Function Description: Removal of foreign material in cavity to access canal orifice

Mode: E  
Power: 1-10

Recommended Service Life  
50 root canals

Water Flow Rate  
mid&high

**EMS E7D**  
**Satelec E7-SD**

**Diamond Coated**

**\$3300**

Function Description: Removal of calcification and filling material in cavity to access canal orifice during re-treatment

Mode: E  
Power: 1-3

Recommended Service Life  
30 root canals

Water Flow Rate  
mid&high

**EMS TC1D**  
**Satelec TC1-SD**

**Diamond Coated**

**\$3300**

Function Description: Removal cervical caries and caries on the occlusal surface

Mode: P  
Power: 1-10

Recommended Service Life  
30 Teeth

Water Flow Rate  
low, mid&high

**EMS E8**  
**Satelec E8-S**

**\$2200**

Function Description: Used to brush the canal walls during re-treatment

Mode: E  
Power: 1-6

Recommended Service Life  
50 root canals

Water Flow Rate  
low, mid&high

**EMS E8D**  
**Satelec E8-SD**

**Diamond Coated**

**\$3300**

Function Description: Used for root canal retrogression, efficient root apical polishing

Mode: E  
Power: 1-6

Recommended Service Life  
50 root canals

Water Flow Rate  
low, mid&high



**Buy 10  
Get 10%  
OFF**  
Mix n Match all  
Scaler Tips



**SCAN ME**  
For more information on  
Scaler Tips

# TOMIDENT **DENTAL SUPPLIES**

## Head Office

96 Elizabeth St.  
Ashfield NSW 2132  
02.9799.1414

[tomident.com](http://tomident.com)

